

GRIP CHIPS™

Grip Chips communicate with compatible i-Limb[®] and i-Digits™ products* via Bluetooth. Users can access automated grips, my grips or sets of favorites when the hand is in close proximity.

The Grip Chips will be delivered preprogrammed with popular grips, as an out of the box ready for use solution for the user. This can be changed to the users needs simply by using the My i-Limb app allowing for the ultimate freedom in personalization.

FEATURES:

- Tapping feature: Double tap on the Grip Chip to access a preprogrammed grip
- Proximity control: Bring the hand in close range to the Grip Chip to access preprogrammed grips
- Ready to use: Preprogrammed with popular grips - no set up required
- Customizable: Use the My i-Limb app to choose the grips available on the Grip Chip



Proximity Control

ÖSSUR TOUCH SOLUTIONS

Össur Touch Solutions are designed to ensure the best possible outcomes for people with upper-limb deficiencies.

* Compatible Devices:
i-Limb Quantum; i-Digits Quantum



Technical Information	
Water Protection:	IP65
Dimensions:	1" x 0,25"
Warranty:	1 year (if sold as part of a package, warranty period of the package will apply.)
Service:	No specific service required. Battery change might be required.
Battery:	CR1616

FOR i-LIMB

Item-No.	Description
TBX50428	PK i-Limb NA Quantum QWD
TBX50429	PK i-Limb NA Quantum WD
TBX50430	PK i-Limb NA Quantum Ti QWD
TBX50431	PK i-Limb NA Quantum Ti WD
TBX50432	PK i-Limb NA Quantum Flexion
TBX50433	PK i-Limb NA Quantum Friction
TBX50434	PK i-Limb NA Quantum TI Flexion
TBX50435	PK i-Limb NA Quantum TI Friction

FOR i-DIGITS

Item-No.	Description
TBX50008	PK1 i-Digits Quantum 1 Digit
TBX50016	PK2 i-Digits Quantum 2 Digit
TBX50024	PK3 i-Digits Quantum 3 Digit
TBX50032	PK4 i-Digits Quantum 4 Digit
TBX50040	PK5 i-Digits Quantum 5 Digit

STANDALONE

Item-No.	Description
PL238300A	Grip Chip kit

FOLLOW ÖSSUR ON



TEL (800) 233-6263
 FAX (800) 831-3160
 WWW.OSSUR.COM

